



Storleksguide gåstolar Miniwalk och Meywalk

	Miniwalk Storlek 1	Miniwalk Storlek 2	Meywalk 2000 Small	Meywalk 2000 Medium	Meywalk 2000 Large	Meywalk MK4 Micro	Meywalk MK4 Mini	Meywalk MK4 Small	Meywalk MK4 Medium	Meywalk MK4 Large
Längd, cm	74	74	85	95	95	77	80	85	86	96
Bredd, cm	60	60	68	68	68	51	56	62	62	66
Bredd insida, cm	39	39	48	48	48	42	44	49	48	53
Sadelhöjd	27-45 *	27-60*	44-62 *	62-86	70-94	25-45	28-52	50-67	64-89	82,5-107
Bålstödshöjd	52-69 *	67-84 *	77-96 *	94-117	108-131	52-65	55-77	74-99	102-135	122-155
Framhjul, mm	150	150	150	200	200	150	150	150	200	200
Bakhjul, mm	250	250	400	400	400	275	275	275	275	275
Vikt, kg	15	16	22	25	26	13	20	25	32	34
Max brukarvikt kg	40	40	100	100	100	60	80	100	130	130