



Storleksguide gåstolar Miniwalk och Meywalk

Mått i cm där inget annat anges	Miniwalk Storlek 1	Miniwalk Storlek 2	Meywalk 2000 Small	Meywalk 2000 Medium	Meywalk 2000 Large	Meywalk 4 Micro	Meywalk 4 Small	Meywalk 4 Medium	Meywalk 4 Large
Längd	74	74	85	95	95	77	85	86	96
Bredd	60	60	68	68	68	51	62	62	66
Bredd insida	39	39	48	48	48	42	49	48	53
Sadelhöjd	27-45	27-60	44-62	62-86	70-94	25-45	44-67	60-84	80-104
Bålstödshöjd	52-69	67-84	77-96	94-117	108-131	52-65	74-99	101-133	120-152
Framhjul, mm	150	150	150	200	200	150	150	200	200
Bakhjul, mm	250	250	400	400	400	275	275	275	275
Vikt, kg	15	16	22	25	26	13	25	32	34
Max brukarvikt kg	40	40	100	100	100	60	100	130	130