

**ECE Group II - III**

15 - 36 kg

33.1 - 79.4 lb

# NAVIGATOR

The world's first foldable rehab car seat



# THOMASHILFEN NAVIGATOR

“ Families who travel a lot want to be flexible and independent. Especially when travelling, a standard rehab car seat is often too unwieldy. The unique THOMASHILFEN **Navigator** as a second seat facilitates the active family everyday life. It can be folded up and taken with you in the twinkling of an eye. Children not only sit safely in the car, but also wherever the environment does not suit them. “



Smart **FOLDING DESIGN** fulfils parents' wishes for small storage space and easy handling on the move.



## CUP HOLDER

Large cup holder (can be used right or left) to keep drinks, snacks and toys ready to hand

## UPHOLSTERY

Cosy upholstery, the breathable cover can be completely removed and washed at 40° C.

## ABDUCTION BLOCK

Soft abduction block adjustable in two positions prevents legs from crossing over

## HEADREST

7-fold height adjustable – easy and quick to adjust



## 5-POINT-BELT

Secure hold with the 5-point positioning belt, practical chest clip prevents twisting of the belts

## LATERAL SIDE SUPPORTS

Additional large lateral side pads optimally support your child



## CARRYING HANDLE

Folded seat can be held comfortably by the carrying handle

## FASTENING BELTS

Integrated fastening belts when the seat is not used in the car but on the road




**Thomashilfen**

In everyday use (examples)



Used as a second seat in the car when the standard rehab car seat is just too bulky or in the other car.



Can also be used on the road in buses and trains if a secure hold while sitting is important.



Sitting safely even on family outings when a well-deserved break is taken.



**You Tube**

Watch the video:  
THOMASHILFEN Navigator  
in everyday use



# THOMASHILFEN Navigator



All details online:  
[www.thomashilfen.com/navigator](http://www.thomashilfen.com/navigator)



## Rehab equipment (selection)



Integrated 5-point positioning belt, practical chest clip prevents twisting of the belts



Soft abduction block adjustable in two positions prevents legs from crossing over



Additional large lateral trunk pads optimally support your child



Fastening belt, if the seat is not used in the car, but on the way

## Accessories (selection)

## 5 upholstery designs available



In the practical transportation bag, the folded seat is comfortably carried over the shoulder or securely stowed away.



GREY



RED



BLUE



GREEN



MAGENTA

Soft, breathable cover for pleasant sitting comfort without annoying sweating. The covers are removable and washable at 40°C.

**You will find complete overview of all accessories on the order form or our website:**  
[www.thomashilfen.com](http://www.thomashilfen.com)

## Dimensions & weights THOMASHILFEN Navigator

Seat depth	32 cm / 12.6"
Seat width front	26 cm / 10.2"
Seat width back	26 cm / 10.2"
Back height (seat)	61 - 71 cm / 24 - 28"
Shoulder height (belt guide)	37 - 47 cm / 14.6 - 18.5"
Folded dimensions (w x h x d)	40 x 43 x 31 cm / 15.7 x 16.9 x 12.2"
Overall dimensions (w x h x d)	68 x 45 x 37 cm / 26.8 x 17.7 x 14.6"
Weight of seat	ca. 5 kg / approx. 11 lb
max. user weight	15 - 36 kg / 33.1 - 79.4 lb

## Basic configuration:

Seat and back with 7-fold height adjustable headrest; complete seat foldable; seat cushion with abduction block (2-fold adjustable); back cushion with side pads; head cushion; additional head cushions (left, right); 5-point positioning belt; chest clip; cupholder (applicable right or left); fastening belt for use outside the car; removable covers, can be machine-washed at 40°C (104°F)

## Navigator

Rückhalteeinrichtung  
für spezielle  
Anforderungen  
Restraint system  
for special needs

S

15 - 36 kg  
33.1 - 79.4 lb



Thomas Hilfen  
GmbH & Co. Medico KG  
27432 Bremervörde  
Germany



Thomas Hilfen für Körperbehinderte GmbH & Co. Medico KG · Walkmühlenstraße 1 · D - 27432 Bremervörde  
 Phone: +49 (0)4761 886-63 or -68 · Fax: +49 (0)4761 886-19 · [www.thomashilfen.com](http://www.thomashilfen.com)