

HOGGI[®]



FLUX 4 walking aid

Spare Parts Catalog

GB/US

When ordering spare parts, please always indicate the following basic data to ensure smooth delivery:

- 1) Serial no. of the product
- 2) Frame width
- 3) Frame colour

Article numbers which have (-xxx) at the end must be completed with a code number for the frame color.

Frame colors	Color code
Grey aluminium (RAL 9007)	(-021)


Contents	Page
Index _____	3
Component Group Overview _____	4
Spare Parts _____	4 - 14
Spare Parts for Accessories _____	15 - 28

Index

Component Group / Spare Parts	Page
Component Group	4
Front frame	5
Rear frame	6
Multifunctional part with revers roll lock left / righth	7 / 8
Stützaufnahme	9
Glider	10
Grab Rail standard	11
Grab Rail universal handle	12
PU-wheels	13
Lower arm support with handles	14
Lower arm support with handles, 3-axis adjustable	15
Hand disk	16
Bumber wheels	17
Drag brake	18
Brake with wheel lock	19
Anti-tip	20
Anti-tip with wheels	21
Seat	22
Hip-/bottom pads	23
Dynamic hip control	24
Carrier	25
Swivel lock	26
Sling seat	27
Sling seat rail	28


Component Group

FLUX Size 4




Item	Article number	Component Group	Page
1-8	3313-0000-021	FLUX Walking aid	
1	3313-1100-021	Front frame	5
2	3313-1200-021	Rear frame	6
3	3313-1800	Multifunctional part with revers roll lock	7 - 8
4	3313-1300-021	Stützaufnahme	9
5	3313-1400	Glider	10
6		Grab Rail	12 - 13
7		Front wheels	14
8		Rear wheels	14
		Accessories	15 - 28

Spare parts





Item	Article number	Article description	Quantity	Unit
1-11	3313-1100E-021	Front frame complete	1	Set.
1-7	3313-1110E-021	Front frame with fasteners	1	Pc.
1	3313-1110-021	Front frame size 4	1	Pc.
2	1062-0412	Blind rivet St/St DIN/ISO 7337 F4x12	4	Pc.
3	1112-0010	Spacer PA black 8x4,3x3	4	Pc.
4	1029-6001	Button head socket screw, zp, Tufflock DIN/ISO 7380 M6x12	8	Pc.
5	1114-0056	Threaded bushing 8xM6x67,5	2	Pc.
6	1114-0057	Threaded bushing 8xM6x70	2	Pc.
7	3313-1160	Hinge bandage D30	4	Pc.
8	1029-0620	Button head socket screw zp DIN/ISO 7380 M6x20	2	Pc.
9	3301-1120	Plug	2	Pc.
10	3301-1130	Saddle bushing	2	Pc.
11	3202-1151	Cover ring	2	Pc.

Spare parts




Item	Article number	Article description	Quantity	Unit
1-12	3313-1200E-021	Rear frame complete grey aluminium	1	Pc.
1-6	3313-1210E-021	Bar, back, with fasteners, grey aluminium	1	Pc.
1	1029-0845	Button head socket screw zp DIN/ISO 7380 M8x45	2	Pc.
2	3201-1270	Plug	2	Pc.
3	3313-1210-021	Bar, back, size 4, grey aluminium	2	Pc.
4	3313-1220	Connector	2	Pc.
5	3313-1240	Hinge bandage oval	2	Pc.
6	3313-1810	Multifunctional part with revers roll lock left	1	Pc.
7	3313-1820	Multifunctional part with revers roll lock right	1	Pc.
8+9	3313-1270E-021	Roll bar	1	Pc.
9	1129-0010	Saddle clamp d 29,8 mm black	2	Pc.
10	3313-1230	Stub shaft	2	Pc.
11	1051-0008	Washer A8,4 zp DIN/ISO 125 A8,4	2	Pc.

Spare parts




Item	Article number	Article description	Quantity	Unit
1-22	3313-1810E	Multifunctional part with revers roll lock left	1	Pc.
1	3313-1251	Revers roll lock bracket	1	Pc.
2	1129-0003	Key ring 24x1,5	1	Pc.
3	3201-7520	Spring housing	1	Pc.
4	1131-0036	Compression spring D-127d 0,8x8,3x29	1	Pc.
5	3313-1817	Plunger block	1	Pc.
6	1039-0001	Screw for plastic, cross rcessed pan head w. flange zp 4x10	1	Pc.
7	1055-0005	Washer A5,3 zp DIN/ISO 9021 A5,3	1	Pc.
8	1042-0006	Hexagon nut self-locking zp DIN/ISO 985 M6	1	Pc.
9	1055-0006	Washer zp A6,4 zp DIN/ISO 9021 A6,4	2	Pc.
10	3313-1253	Shaft axle	1	Pc.
11	1029-0630	Button head socket screw zp DIN/ISO 7380 M6x30	1	Pc.
12	1080-0010	Retaining ring for shaft DIN/ISO 471 10x1	1	Pc.
13	1059-0001	Shim ring zp DIN/ISO 988 10x16x1,2	1	Pc.
14	3301-1830	Freewheel left	1	Pc.
15	3301-1812	Freewheel axle	1	Pc.
16	3301-1814	Spacer bushing	1	Pc.
17	1029-6001	Button head socket screw, zp, Tufflock DIN/ISO 7380 M6x12	2	Pc.
18	1024-0616	Hexagon socket countersunk head screw zp DIN/ISO 7991 M6x16	1	Pc.
19	3301-1816	Spring hanger	1	Pc.
20	1133-0012	Leg spring T-18854L 1,6x14,4x3,00	1	Pc.
21	3301-1811-001	Base plate left	1	Pc.
22	3313-1815	Axis	1	Pc.

Spare parts




Item	Article number	Article description	Quantity	Unit
1-22	3313-1820E	Multifunctional part with revers roll lock rigth	1	Pc.
1	3313-1251	Revers roll lock bracket	1	Pc.
2	1129-0003	Key ring 24x1,5	1	Pc.
3	3201-7520	Spring housing	1	Pc.
4	1131-0036	Compression spring D-127d 0,8x8,3x29	1	Pc.
5	3313-1817	Plunger block	1	Pc.
6	1039-0001	Screw for plastic, cross rcessed pan head w. flange zp 4x10	1	Pc.
7	1055-0005	Washer A5,3 zp DIN/ISO 9021 A5,3	1	Pc.
8	1042-0006	Hexagon nut self-locking zp DIN/ISO 985 M6	1	Pc.
9	1055-0006	Washer zp A6,4 zp DIN/ISO 9021 A6,4	2	Pc.
10	3313-1253	Shaft axle	1	Pc.
11	1029-0630	Linsenschraube M6x30	1	Pc.
12	1080-0010	Retaining ring for shaft DIN/ISO 471 10x1	1	Pc.
13	1059-0001	Shim ring zp DIN/ISO 988 10x16x1,2	1	Pc.
14	3301-1840	Freewheel right	1	Pc.
15	3301-1812	Freewheel axle	1	Pc.
16	3301-1814	Spacer bushing	1	Pc.
17	1029-6001	Button head socket screw, zp,Tufflock DIN/ISO 7380 M6x12	2	Pc.
18	1024-0616	Hexagon socket countersunk head screw zp DIN/ISO 7991 M6x16	1	Pc.
19	3301-1816	Spring hanger	1	Pc.
20	1133-0013	Leg spring T-18854L 1,6x14,4x3,00	1	Pc.
21	3301-1821-001	Base plate right	1	Pc.
22	3313-1815	Axis	1	Pc.

Spare parts




Item	Article number	Article description	Quantity	Unit
1-7	3313-1300E-021	Stützaufnahmen	2	Pc.
1	3301-1240	Clamp bushing	1	Pc.
2	1129-0007	Saddle clamp with quick-release D 31,8 left	1	Pc.
	1129-0006	Saddle clamp with quick-release D 31,8 right	1	Pc.
3	3213-1400E	Glider complete	1	Pc.
4	3313-1310-021	Stützaufnahme	1	Pc.
5	3313-1331	Glider stopper	1	Pc.
6	3313-1320	Support tube bearing	1	Pc.
7	3313-1232E	Hakenfang with fasteners	1	Pc.

Spare parts




Item	Article number	Article description	Quantity	Unit
1-5	3213-1400E	Glider complete	1	Pc.
1	1029-0612	Button head socket screw zp DIN/ISO 7380 M6x12	4	Pc.
2	1114-0028	Threaded bushing 8xM6x38	2	Pc.
3	1133-0016	Leg spring T-18805L 1,6x12, 1x3,25	1	Pc.
4	3313-1411	Glider	1	Pc.
5	3313-1412	Safety hook	1	Pc.

Spare parts



Item	Article number	Article description	Quantity	Unit
1-10	3313-1500E	Grab rail standard, size 4 (retrofit)	1	Pc.
1-6	3313-1510E-006	Bracket left complete, anodised bright	1	St.
2-7	3313-1520E-006	Bracket right complete, anodised bright	1	St.
2	1134-0006	Snap button B 141 7,7x9,4	1	St.
3	1141-0015	Ribbed insert GL 25x1-2.5	1	St.
4	3301-1240	Clamp bushing	1	St.
5	1129-0008	Saddle clamp with bolt D 31,8	1	St.
6	3301-1530	Flange bushing	1	St.
8-10	3313-1540	Handle standard	1	St.
8	3313-1541-006	Tube bend handle standard size 4 anodised bright	1	St.
9	1062-0408	Blind rivet St/St DIN/ISO 7337 F4x8	2	St.
10	1143-0025	Handle grip TPE black d25 x 120	2	St.

Spare parts



Item	Article number	Article description	Quantity	Unit
1-14	3313-5100E	Grab rail with universal handle, size 4 (retrofit)	1	Pc.
1-6	3313-1510E-006	Bracket left complete, anodised bright	1	Pc.
2-7	3313-1520E-006	Bracket right complete, anodised bright	1	Pc.
2	1134-0006	Snap button B 141 7,7x9,4	1	Pc.
3	1141-0015	Ribbed insert GL 25x1-2.5	1	Pc.
4	3301-1240	Clamp bushing	1	Pc.
5	1129-0008	Saddle clamp with bolt D 31,8	1	Pc.
6	3301-1530	Flange bushing	1	Pc.
8	3303-5110-006	Tube bend handle universal anodised bright	1	Pc.
9-13	3301-5120	Universal handle complete	2	Pc.
9	3301-5130-006	Handle universal anodised bright	1	Pc.
10	1089-0002	Wedge M8x21	1	Pc.
11	1021-0885	Hexagon socket head cap screw zp DIN/ISO 912 M8x85	1	Pc.
12	1051-0008	Washer A8,4 zp DIN/ISO 125 A8,4	1	Pc.
13	1070-0308	Spring pin DIN/ISO 7346 3x8	1	Pc.
14	1143-0005	Handle grip TPE black d25 x 120	1	Pc.

Spare parts



Item	Article number	Article description	Quantity	Unit
1	1281-0007	Castor fork LG-K-200-12 black 200x60	2	Pc.
2	1205-0875	Hexagon head screw M8x75 with nut	2	Pc.
3	3301-1150	Steering pivot pin	2	Pc.
4	1051-0012	Washer A13 zp DIN/ISO 125 A 13	2	Pc.
5	1056-0013	Toothed lock washer A13 zp DIN/ISO 6797A A12,5	2	Pc.
6	1221-2048	PU-Wheel R200/48 (chrome rim)	2	Pc.
7	1221-1950	PU-wheel 190x50	2	Pc.
8	1221-3061	PU-Wheel R300/60-6-LSSI-BS-SV (chrome) with break disk	2	Pc.
9	1221-2961	PU-wheel 290x60 with break disk	2	Pc.

Reservdelslista FLUX stl 4

Art.nr	Artikel
1021-0625	Sexkanthålskruv helgängad M6x25
1021-0635	Sexkanthålskruv DIN 912 M6x35
1021-0675	Sexkanthålskruv DIN/ISO 912 M6x75
1021-0820	Sexkanthålskruv DIN/ISO 912 M8x20
1021-0885	Sexkanthålskruv DIN 912 M8x85
1022-0625	Sexkanthålskruv DIN/ISO 7984 M6x25
1024-0616	Sexkantskruv DIN/ISO 7991 M6x16
1029-0612	Sexkantshålskruv med rundat huvud zp DIN/ISO 7380 M6x12
1029-0620	Sexkantshålskruv med rundat huvud zp DIN/ISO 7380 M6x20
1029-0630	Sexkantshålskruv med rundat huvud DIN/ISO 7380 M6x30
1029-0645	Sexkantshålskruv med rundat huvud zp DIN/ISO 7380 M6x45
1029-0650	Sexkantshålskruv med rundat huvud zp DIN/ISO 7380 M6x50
1029-0655	Sexkanthålskruv med rundat huvud DIN/ISO 7380 M6x55
1029-0845	Sexkanthålskruv med rundat huvud DIN/ISO 7380 M8x45
1029-6001	Skruv med rundat huvud,Tufflock DIN/ISO 7380 M6x12
1031-0830	Sexkantshålskruv med rundat huvud zp DIN/ISO 913 M8x30
1039-0001	Skruv för plast 4x10
1042-0006	Sexkantslåsmutter DIN/ISO 985 M6
1042-0008	Sexkanthålsmutter DIN/ISO 985 M8
1051-0006	Bricka A6,4 zp DIN/ISO 125 A6,4
1051-0008	Bricka A8,4 zp DIN/ISO 125 A8,4
1051-0012	Bricka A13 zp DIN/ISO 125 A 13
1055-0005	Bricka A5,3 zp DIN/ISO 9021
1055-0006	Bricka zp A6,4 zp DIN/ISO 9021 A6,4
1055-0008	Bricka zp A8,4 zp DIN/ISO 9021 A8,4
1056-0013	Bricka A13 zp DIN/ISO 6797A A12,5
1059-0001	Passbricka DIN/ISO 988 10x16x1,2
1062-0408	Blindnit St/St DIN/ISO 7337 F4x8
1062-0412	Blindnit DIN/ISO 7337 F4x12
1070-0308	Spår sprint DIN 7346 3x8
1080-0010	Spårring för axel DIN/ISO 471 10x1
1089-0002	Metall kil M8
1089-0004	Låsband 525 MXR
1099-0011	Skruv i vinkel, 160° angle
1099-0015	Nippel
1111-0001	Distans kloss, svart PA 25x6.3x8
1112-0010	Distansbricka, svart PA 8x4,3x3

1113-0047	Brica, svart, 24x8,4x4
1114-0028	Bussning 8xM6x38
1114-0056	Bussning 8xM6x67,5
1114-0057	Bussning 8xM6x71
1116-0003	Rörklämma, svart PA D30
1116-0005	Rörklämma, svart PA klädd insida D25
1127-0009	Låssprint M10
1127-0010	Låssprint ZS-10x1-OA-KS d6xM10xS10x31 (w/o catch)
1129-0003	Nyckelring 20x1,5
1129-0006	Rörklämma med snabbås D 31,8 höger
1129-0007	Rörklämma med snabb lås D 31,8 vänster
1129-0008	Rörklämma med lås D 31,8
1129-0010	Rörklämma 29,8 mm black
1131-0034	Fjäder D-180Q-03 1,2x11x46,7
1131-0036	Tryckfjäder D-127d 0,8x8,3x29
1131-0045	Fjäder D-207J-90 1,5x18,5x104
1131-0046	Fjäder D-251 2,0x18x98
1133-0012	VridfjädrarT-18855L 1,6x14,4x3,25
1133-0013	Vridfjädrar T-18855R 1,6x14,4x3,25
1133-0016	Vridfjäder T-18805L 1,6x12, 1x3,25
1134-0006	Knapplås med fjäder B 141 7,7x9,4
1134-0007	Knapplås D-150
1141-0006	Ändavslutning SFL 20x0,8-2,5
1141-0009	Ändavslutning SFL 25x1-2,5
1141-0015	Ändavslutning GL 25x1-2.5
1142-0002	Huv utvändig, svart K 20 PE
1143-0005	Handtag, svart TPE d25 x 120
1143-0006	Mjukt skumrör, svart TP/K
1145-0001-015	Connector T-90° 25/25 with bore 4,2
1145-0002-015	Connector T-90° 25/25
1150-0001	O-ring 18x2
1150-0002	O-ring NBR 70 DIN/ISO 3771 10x2
1204-0005	Bromsdisk i aluminium
1205-0875	Skruv M8x75
1209-0002	Handbroms 80° manuellt lås 25mm
1209-0019	Kabel 1300 / 1500
1221-1950	PU-hjul 190x50
1221-2048	PU-hjul R200/48 (chrome rim)
1221-2961	PU hjul 290x60 med bromsdisk
1221-3061	PU-hjul R300/60-HOG-BS-SV (chrome) with break disk
1249-0001	Avbärrarhjul AR75/15,5/9005 black
1249-0002	Avbärrarhjul 100 mm black
1249-0008	Hjul på tippstöd
1281-0007	Gaffel LG-K-200-12 schw 200x60
1511-2535	Strap 25x350
1512-0006	Strap med fäste
1512-0007	Strap med fäste
1530-0004	Fäste
3201-1270	plug
3201-7520	Fjäderhylsa

3202-1151	Bussning
3202-7100	Riktningsspärr för ledade länkhjul fram, set
3202-7110E	Riktningsspärrhylsa med sexkantsskruv
3301-1120	Gaffelhylsa
3301-1130	Sadel bussning
3301-1150	Bultaxel
3301-1240	Klammerbussning
3301-1530	Kopplingsbeslag
3301-1811-001	Basplatta vänster sida
3301-1812	Frihjulsaxel
3301-1814	Distansbussning
3301-1816	Fjäderhylsa
3301-1821-001	Basplatta höger
3301-5120	Handtag komplett sats för universalhandtag
3301-5210	Terminal holder
3301-5211	Dubbelbeslag 60
3301-5220-006	Handtag, svart
3301-5221-006	Handtag, anodized
3301-5321-006	Handhold
3301-5360	Distans
3301-5410	Sadelbussning
3301-5500	Släpbroms, set
3301-5510	Släpbromsplatta
3301-5520	Gängad bussning
3301-5810	Sittplatta
3301-5910	Klammerbeslag komplett
3301-5913	Axel med genomgående hål
3301-5915	Klammerbeslag 45
3301-5920	Support arm pelotte
3301-5921	Rör, endast
3301-5933	Distans
3301-5950	Centerdyna
3301-5960	Höftdyna
3302-5240	Dyna
3303-5110-006	Bakre båge till universalhandtagsbågen stl 3
3303-6110	Pakethållare stl 3 endast rör
3313-1100E-021	Främre ram, komplett
3313-1110-021	Främre ram, endast rör
3313-1110E-021	Främre ram
3313-1160	Plugg
3313-1200E-021	Bakre ram, komplett
3313-1210-021	Sidorör
3313-1210E-021	Bakre ram
3313-1220	Svart plast fäste (till tippstöd)
3313-1230	Fäste
3313-1240	hinge bandage oval
3313-1251	Fäste till backspärr
3313-1253	Axel
3313-1270E-021	Roll bar
3313-1300E	Sidorör, komplett

3313-1310-021	Sidorör
3313-1320	support tube bearing
3313-1331	Säkerhetsstopp
3313-1332E	Fäste
3313-1400E	Fällningsmekanism, komplett
3313-1411	Fällningsmekanism
3313-1412	Säkerhetsspärr
3313-1510E-006	Höjdjusteringsrör, komplett, vänster
3313-1520E-006	Höjdjusteringsrör, komplett, höger
3313-1541-006	Båge stl 4
3313-1810	Multifunktionellt fäste, vänster
3313-1815	Axel
3313-1817	Plunger block
3313-1820	Multifunktionellt fäste, höger
3313-1830	Frihjul, vänster
3313-1840	Frihjul, höger
3313-5100	Universalhandtag
7000-0022	Gummiband med krokar